


natural-fill

FOR NATURAL BEAUTY

Cannula can be performed through a single access point, and the extended length of a cannula allows for the selection of the safest incision area. Also it can be bent around tissue structures and push away blood vessels and nerve, where a needle would simply pass through. Thus, cannula can reduce pain and shorten recovery time, making the procedure more comfortable for the patient. Furthermore, blunt cannula helps to reduce the risks of intravenous arterial embolism, enhancing patient safety.

- Less downtime • Filling the solution
- Liposuction • Safety, do not injure the vein
- Filling the several places with one cannula



CANNULA SPEC

Anatomical Region	Cannula Recommendation		
Periorbital	25, 27, 30G : 25 ~ 40mm	17G	70 / 90mm
Forehead	25, 27G : 38 ~ 50mm	18G	50 / 60 / 70 / 90mm
Brow	25, 27G : 25 ~ 50mm	19G	50 / 60 / 70 / 90mm
Temporal, Middle third of the face, Lower third of the face	25, 27G : 38 ~ 50mm	20G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
Lips	25 - 30G : 25 ~ 40mm	21G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
Chin, Neck, Decollete, Hands	25 - 27G : 38 ~ 50mm	22G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
		23G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
		24G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
		25G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
		26G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
		27G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
		29G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm
		30G	25 / 30 / 35 / 40 / 50 / 60 / 70 / 90mm