



# natural-meso

FOR NATURAL BEAUTY



The thin wall needle enables less extrusion force and results in smoother injection. Meso therapy needle is effectively used for various ways such as injecting hyaluronic acid, local anesthetics, mesotherapy and placenta as well as high viscous solution.

27G	4mm
	13mm
30G	4mm
	6mm
	8mm
	10mm
	13mm
31G	4mm
	6mm
	8mm
	10mm
	13mm
32G	4mm
	6mm
	8mm
	10mm
	13mm
33G	4mm
	6mm
	8mm
	13mm
34G	4mm
	6mm
	8mm
	13mm